



Cleaning Tool


For automated cleaning



For manual cleaning



For purging



Cleaning Tube

	order code PD-CTTU-606501 Cleaning Tube 6.0/L650 <small>Supplied with Cleaning Holder PD-CTH</small>		order code PD-CTTU-600504 Cleaning Tube 6.0/L50 <small>Supplied with Cleaning Holder PD-CTH</small>		order code PD-CTTU-150252 Cleaning Tube 1.5/L25 <small>Supplied with Cleaning Adaptor IN PD-CTAD-IN</small>
--	--	---	---	---	---

Spray Adaptor for Slim Motor Handpieces/Attachments (P200-EZ-SM)

P300 Attachment

- Slim Attachment 300
- Standard Attachment 300

order code
PD-EZ-M
EZ Spray Adaptor

- PANA Spray Plus is not included

order code
Z182600
PANA SPRAY Plus



CREATE IT.

NAKANISHI INC.
700 Shimohinata, Kanuma Tochigi 322-8666, Japan
TEL : +81 289 64 3380 FAX : +81 289 62 5636

NSK America Corp.
1800 Global Parkway, Hoffman Estates, IL 60192, USA
TEL : +1 888 675 1675 FAX : +1 800 838 9328

Primado2
Slim Attachment 300




Slim Attachment 300 80,000min⁻¹

Slim Attachment 300 is designed with new features to deliver improved visibility, smooth operation, safety, and stable rotation performance.

Variable exposure

· Adjustable exposure lengths of 7/9/11/13/15/17mm.




Non-slip surface

· Improved grip with surface etching.



Accessory

· Two types of Guards can slide on to the attachment



Quick Joint

· Easy attachment/detachment with Quick Joint
· 16 different handpiece combinations are possible with the straight/angled hubs and variety of tubes.



Smooth connection

· Quick Joint allows for an easy attachment-motor connection.



Tapered Tip

· Ensures better visibility and access to surgical anatomy.



Compatible with Primado 2 Control Unit and Slim Motor



order code
P200-CU-120
Primado2 Control Unit for 120V



order code
FC-73/FC-74
Primado2 Single Foot Control
Primado2 Multi Foot Control




order code
P200-SMH-HS
Primado2 Slim Motor/HS
· Hand Control Motor
· Speed : Max. 80,000 min⁻¹

order code
P200-SMH
Primado2 Slim Motor
· Speed : Max. 80,000 min⁻¹

Operability

Straight attachment added to the product lineup.



Reduce hand fatigue with the Straight Hub.



By using a Straight Hub, surgeons can easily access a narrow space by simple rotation.



Non-slip surface etching improves grip.

Accessibility


Better accessibility with the 6 exposure lengths adjusting mechanism. The length can be adjusted with the Accessory attached.




Visibility

Visibility is improved with the tapered tip and variable exposure lengths.

exposure length 11mm



exposure length 17mm



Lineup

Slim Attachment Hub




order code
P300-1AHS
Slim Attachment Hub/Straight



order code
P300-1AHA
Slim Attachment Hub/Angled

Slim Tube



order code
P300-1T110
Slim Tube 110



order code
P300-1T110-C
Slim Tube 110/C



order code
P300-1T110-PB
Slim Tube 110/PB



order code
P300-1T130
Slim Tube 130




order code
P300-1T130-C
Slim Tube 130/C



order code
P300-1T130-PB
Slim Tube 130/PB



order code
P300-1T170
Slim Tube 170



order code
P300-1T240
Slim Tube 240

Accessory



order code
P300-1T-B



order code
P300-1T-H

*Cannot be used with an irrigation nozzle.

Irrigation Nozzle

















order code
P200-IN-1A155
For P300-1T110 / P300-1T110-C / P300-1T110-PB

order code
P200-IN-1A175
For P300-1T130 / P300-1T130-C / P300-1T130-PB



Cutting Accessories for Slim Attachment 300

Cutting Accessories								Head size ● (mm) Actual size	Accessory	Slim Tube	Slim Attachment Hub
R Fluted Burs	R Diamond Burs	R Coarse Diamond Burs	MH Diamond Burs	MH Coarse Diamond Burs	RT Diamond Burs	RT Coarse Diamond Burs					
order code	order code	order code	order code	order code	order code	order code	order code				
—	PDS-1RD110-20	PDS-1RCD110-20	PDS-1MHD110-20	—	—	—	—	● 2.0	 order code P300-1T-B	110 mm	
PDS-1RF110-30	PDS-1RD110-30	PDS-1RCD110-30	PDS-1MHD110-30	PDS-1MHCD110-30	—	—	—	● 3.0	 order code P300-1T-H	 order code P300-1T110 Slim Tube 110	
PDS-1RF110-40	PDS-1RD110-40	PDS-1RCD110-40	PDS-1MHD110-40	PDS-1MHCD110-40	—	—	—	● 4.0	—	 order code P300-1T110-C Slim Tube 110/C	
—	—	—	—	—	PDS-1RTD110-45	PDS-1RTCD110-45	—	● 4.5	—	 order code P300-1T110-PB Slim Tube 110/PB	
—	—	PDS-1RCD110-50	—	—	—	—	—	● 5.0	—	—	
—	PDS-1RD130-20	PDS-1RCD130-20	PDS-1MHD130-20	—	—	—	—	● 2.0	 order code P300-1T-B	130 mm	
PDS-1RF130-30	PDS-1RD130-30	PDS-1RCD130-30	PDS-1MHD130-30	PDS-1MHCD130-30	—	—	—	● 3.0	 order code P300-1T-H	 order code P300-1T130 Slim Tube 130	
PDS-1RF130-40	PDS-1RD130-40	PDS-1RCD130-40	PDS-1MHD130-40	PDS-1MHCD130-40	—	—	—	● 4.0	—	 order code P300-1T130-C Slim Tube 130/C	
—	—	—	—	—	PDS-1RTD130-45	PDS-1RTCD130-45	—	● 4.5	—	 order code P300-1T130-PB Slim Tube 130/PB	
—	—	PDS-1RCD130-50	—	—	—	—	—	● 5.0	—	—	
PDS-1RF170-35	PDS-1RD170-35	PDS-1RCD170-35	—	—	—	—	—	● 3.5	—	170 mm	
PDS-1RF170-40	PDS-1RD170-40	PDS-1RCD170-40	—	—	—	—	—	● 4.0	—	 order code P300-1T170 Slim Tube 170	
—	PDS-1RD170-50	PDS-1RCD170-50	—	—	—	—	—	● 5.0	—	—	
—	PDS-1RD240-30	PDS-1RCD240-30	—	—	—	—	—	● 3.0	order code P300-1T-B order code P300-1T-H	240mm	
PDS-1RF240-40	PDS-1RD240-40	PDS-1RCD240-40	—	—	—	—	—	● 4.0	—	—	
—	PDS-1RD240-45	PDS-1RCD240-45	—	—	—	—	—	● 4.5	—	 order code P300-1T240 Slim Tube 240	
—	PDS-1RD240-50	PDS-1RCD240-50	—	—	—	—	—	● 5.0	—	—	